

THE PRECISION SCIENTIFIC CO (CBE)



great

# Combo Offer

Wee Series

TFT color capacitive touch screen  
(5 inches, 800x480 pixels)



Gradient PCR (WEE32)  
Electrophoresis (LA666)  
Power Supply (LA690)

HIMEDIA

For Life is Precious



WiFi Enabled  
Ultra light Sleek  
Compact 3kg

Rs. 1,43,050 /-

( GST-Extra )

Capacity	Tube	Temperature range	MAX. Heating Rate	MAX. Cooling Rate	Uniformity	Accuracy	Display Resolution	Temperature Control	Ramping Rate Adjustable
32X0.2ml (4X8 LAYOUT)	0.2ml single tube, 8 strip	4~100° C	5° C/s	4° C/s	±0.2° C	±0.2° C	0.1° C	Block/ tube	0.1° C-0.5° C

## Electrophoresis Power Supply

LA 690  
4 Terminal



### Accessories

Connecting Cable : 1 Set (Black & Red)  
Output Range : 10 to 300 VDC

### Power Supplies

No of Outputs : 4  
Compatible with : all Electrophoresis range

## Hi-Gel Run1014

LA 666



### Accessories

Gel Casting Tray : 1 No  
Well Comb : 2 Nos. (10 well & 14 well)  
Detachable Electrodes (Anode, Cathode)  
Connecting Cable : 1 Set (Black & Red)

### Specification

Tank Size (LxWxH) in cm: 22x15x9.5  
Gel Tray Size (in cm) : 9x11  
Sample Capacity (nos) : 24  
Max. Buffer Volume (in ml) : 600  
Comb Configurations (teeth) : 10 & 14  
Compatible Power Supply : La690



THE PRECISION SCIENTIFIC CO., ( CBE )

0422-2235970 , 2237280 , 4506280 , 9786022732

info@pscbe.com