



# Product Catalogue

## MICROSCOPES & ACCESSORIES

#magnusmicroscopes



### ABOUT US

Driven and committed to developing and manufacturing products for today's scientific and medical community. We strive to continuously evolve to stay abreast with the speed of change in the technology world.

# MLXi Plus

## Biological Microscope



- Plan infinity corrected optics
- Multilayer coated optics
- Ergonomic sidentopf head
- FN.20 wide field eyepiece
- Inward facing quadruple revolving nosepiece
- Available in LED, halogen & battery back-up variants
- Anti-fungal optics
- Hard-coated stage with double slide holder
- Rackless stage - optional

Model	Head	Objectives	Eyepiece	Stage	Condenser	Illumination
MLXi Plus	30° Binocular, Sidentopf	4x, 10x, 40x, 100x (oil) Plan	WH 10x (FN. 20)	Hard-coated, low drive mechanical	1.25 Abbe, Rack and Pinion, Iris	6v 20w Halogen
MLXi Plus LED	30° Binocular, Sidentopf	4x, 10x, 40x, 100x (oil) Plan	WH 10x (FN. 20)	Hard-coated, low drive mechanical	1.25 Abbe, Rack and Pinion, Iris	1w & 3w LED
MLXi Plus LED Freedom	30° Binocular, Sidentopf	4x, 10x, 40x, 100x (oil) Plan	WH 10x (FN. 20)	Hard-coated, low drive mechanical	1.25 Abbe, Rack and Pinion, Iris	1w LED with battery back-up

Magnus